


VERIFY Tagged Antigen™

Validation Data

Potential Applications

1. Protein purifications
2. Protein function studies
3. Antibody validations in Western blot, ELISA, immunoprecipitation, etc.
4. Assay standards
5. Protein-protein interaction
6. Prepare reverse phase protein arrays

Protein Production in HEK293T Cells



HEK293T cells were transfected with TrueORF INPP5K clone (RC202190). Cells were lysed 48hrs later for protein purification using M2 beads (anti-DDK). Final purified protein is > 80% pure in Comarssie Blue staining.

Study Protein Functions

Protein purification using anti-tag antibody and *in vitro* function assays

Vol 456 | 20 November 2008 | doi:10.1038/nature07470

nature

ARTICLES

Identification of Holliday junction resolvases from humans and yeast

Stephen C. Y. Ip^{1*}, Ulrich Rass^{1*}, Miguel G. Blanco^{1*}, Helen R. Flynn², J. Mark Skehel² & Stephen C. West¹

- For over 20 years, little was known about eukaryotic Holliday junction resolvase
- GEN1 TrueORF clone with FLAG tag (RC221451) was transfected into HEK293T cells
- GEN1 protein was purified using anti-FLAG M2 affinity column
- Purified GEN1 protein resolved Holliday junction X0 efficiently *in vitro*
- Both GEN1 and Yen1 were identified as resolvases for eukaryotic cells






Figure 3 | Resolution of Holliday junctions by recombinant Yen1 and GEN1. **a**, ³²P-labelled Holliday junction X0 was incubated with cell-free extracts from yeast overexpressing Flag-tagged Yen1, or a catalytically inactive Yen1 (E193A/E195A) mutant (Yen1^{mut}), and the products were analysed by neutral PAGE. Control: extracts from cells transformed with empty expression vector. Yen1-Flag and Yen1^{mut}-Flag were detected by western blotting using anti-Flag antibody. **b**, As **a**, but affinity purified GEN1-Flag, or ResA, was used. **c**, **d**, Holliday junction X0, 5'-³²P-end-labelled in strand 1 or 3, was treated as in **a** and **b**, and products were analysed by denaturing PAGE. Asterisks indicate the strand with the radioactive label. **e**, The cleavage sites. **f**, Wild-type and mutant derivatives of GEN1-Flag were assayed using Holliday junction X0. Lane 1, control lane; lane 2, GEN1(D30A); lane 3, GEN1(E134A); lane 4, GEN1(E136A); lane 5, GEN1(D157A); lane 6, GEN1(E134A/E136A); lane 7, wild-type GEN1. In the bottom panel, GEN1-Flag proteins were detected by western blotting.

Antibody Validation in WB


Over-expression cell lysate of human STAT3 (NM_139276) was used to test 3 commercial antibodies. Antibody A shows strong antigen binding. Antibody B did not react with STAT3 at all, while Antibody C shows weak binding.



Over-expression cell lysate for human LRRK2 (NM_198578) was used to validate two anti-LRRK2 antibodies in two rounds of Western blot experiments. The left panel shows positive reaction with a commercial monoclonal antibody. The right panel shows a negative reaction with a polyclonal antibody generated against a protein peptide.



Antibody Validation – P53 Ab




Five commercial antibodies against human P53 were evaluated in Western blot experiments with P53 over-expression cell lysate.

P53 protein level in cell lysate was pre-determined using a purified GST-Myc-DDK standard. Lysate was serially diluted before SDS-PAGE and immunoblotting.

Antibodies A-C are rabbit antibodies. D and E are mouse antibodies.

Antibody quality and star rating is based on P53 protein detection level.


Determining Lot to Lot Consistency



Three different lots of a commercial polyclonal antibody against human P53 were evaluated in Western blot experiments with P53 over-expression cell lysate.

One can see certain lot to lot variation in amount of antigen detection, especially between lot A and C.

Endogenous Protein Expression Level




Estimate the endogenous p53 level in different carcinoma cell line by VERIFY antigen standard. 10 ug whole cell lysates from 3 different cancer cell lines were tested against titrated p53 over-expression lysate. Using over-expression lysate as standard, it is determined endogenous p53 expression level in these 3 cell lines were 3.8 ng (A431), 0.4 ng (HT-29), and 1.3 ng (HEKT293) respectively.

Reverse Phase Lysate Arrays

Figure shows a preliminary 4,000 protein array results when spotted lysate arrays are detected using anti-FLAG antibody and anti-P53 antibody each labeled with different fluorescent dyes (Alexa-555 and Alexa-647).

Left two panels are individual antibody staining, while right panel shows combined image to illustrate differential binding.

Protein arrays can be used in detection of autoimmune antibodies, protein-protein interaction, and antibody decoding.




Anti-FLAG

Anti-p53

Merged

Over-expression Cell Lysate QC

Figures show typical Western blot images when over-expression cell lysates were tested with anti-tag antibodies. Figure 4 shows a comparison using either anti-myc antibody (TA100010) or anti-DDK antibody (TA100011). Figures 5 and 6 were done with anti-myc antibody. Most proteins are detected as a single band. Some proteins are detected as two or multiple bands due to various reasons, such as post-translation modification, protein degradation, etc.



* Molecular weight is smaller than 10 kDa.

Figure 6

Over-expression Cell Lysate QC




Figure 7




Figure 8